

LaserToF LT2 Plus

## BENCHTOP MALDI-ToF MS

SAI's LaserToF LT2 Plus is a completely automated, compact and versatile matrix assisted laser desorption ionisation time of flight (MALDI-TOF) mass analyser. The instrument is designed for the rapid and accurate analysis of a wide range of biomolecules such as peptides, proteins and oligonucleotides. The LT2 Plus is also suitable for the analysis of industrial polymers as used in the petroleum and chemical additive industries.

# **Key Features**

- Linear MALDI-TOF
- 96 Well 1/4 Microtitre Plate
- Pulsed Extraction
- Simple Operation
- Sample Viewing Camera
- Attractive Pricing
- Small Footprint



# **Mass Spectrometer**

- Advanced acquisition and control electronics
- High stability and easy calibration for optimum precision for a wide range of measurements
- Fast electron multiplier detector system as standard for linear operation
- High gain, large dynamic and linear range to improve the system performance
- Detector has long lifetime and low maintenance



## **Data System**

#### Hardware

Windows PC including:

- 2 Gbytes of RAM
- · High capacity hard disk drive
- CD-ROM drive
- 19" flat screen colour monitor.
- Network connection

### Software

Full control of all MS parameters,

- Data acquisition and data reduction program
- Automated sequence programming for batch operation
- Database search of peptide map data
- The system comes network ready

# **Specifications**

• **Resolution** >2,500 (FWHM)

Mass Accuracy

### Laser

Type Pulsed UV 337nm(nitrogen)

## **Vacuum System**

Pumping is provided by a high efficiency air cooled turbomolecular pump coupled to an oil free diaphragm backing pump that is built into the same footprint.

#### **Services**

Power: 220-250V, 4.5 A, single phase


or 110V, 9A single phase

Environment: Temperature range from 15-27°C Dimensions: 83 x 52 x 130 cm (W x D x H)

Weight: Approx. 130 Kilos

Services: Electricity only. No water required

## Layout



Scientific Analysis Instruments Ltd Hadfield House, Hadfield Street Manchester M16 9FE

Tel: + 44 (0) 161 874 2460 Fax: + 44 (0) 161 874 2461 E-mail:sales@saiman.co.uk http://www.saiman.co.uk **Local Representative** 

SAI has a policy of continual improvements to our products. We believe the information in this document is correct. However, it is subject to change without notice. SAI cannot be held responsible for any errors or omissions in this document.